

ENTFirst

Ear, Nose and Throat Treatment Unit

FU100

Japanese technology and spirit cultivated for over 60 years, distilled into one exceptionally reliable unit.

Made in Japan

Easy to Use

Space Saving

Reliable





Easy to Use

All-in-one unit with casters. Lamp and pole are included as standard. Microscope (option) can be mounted on the lamp pole.

Reliable

Highly reliable mechanism

Space Saving

Can be installed on just 65 cm x 52 cm floor space.

Reliable

Clean, convenient spray (Autoclavable, nozzle detachable from the spray main body with 90 degree rotation.)



Ear, Nose and Throat Treatment Unit

FU100

FU100 fulfils the needs of all medical staff and patients. FU100 is manufactured with pride of "Made in Japan" and our confidence. The simple and smart design puts patients at ease. The outstanding operation and quality assure medical staff satisfaction and safety.

Easy to Use

Simple and direct operation of main functions using three buttons

Space Saving

Slide-in/pull-out lid on surgical instrument storage tray

Space Saving

Retractable hoses



Specifications

Spray device	Cotton-covered silicone tube (antibacterial) / free-joint type (Straight: 2 Curved: 1)	3
Suction device	Oil-free suction motor	1
	Cotton-covered silicone tube (antibacterial) with pulley (Option: Exterior loop type)	1
	Hard glass secretion bottle (3-liter)	1
Ventilation air device	Push-button type	1
Air supply	Oil-free compressor	1
	Stainless steel tank	1
	Automatic drain system	1
Laryngealmirror heater	Hot air blower (with thermostat and temperature fuse)	1
Lamp*	Heat-absorbing glass, double arm type Pole diameter 49 mm	
Unit Material	Body (steel/melamine coating) Top panel (artificial marble)	
Dimensions	650 (W) x 530 (D) x 800 (H) mm (Max 2,450 (H) mm when lamp is attached)	
Weight	94 kg (When lamp is attached)	
Rated voltage	AC 230 V ±10 % 50 Hz	
Consumption current	7 A	

* Lamp can be removed. Specify when ordering.

ENTFirst

DAIICHI MEDICAL CO., LTD.

2-27-16, Hongo, Bunkyo-ku, Tokyo 113-0033, Japan Tel: 81-3-3814-0111 Fax: 81-3-3814-0135 Email: first-ex@first-med.co.jp

<http://www.first-med.co.jp>